



DRY BATH INCUBATOR BDIB-117

DRY BATH INCUBATOR BDIB-117

DRY BATH INCUBATOR

Round metal blocks and more compact structure, suitable for most centrifuge tubes. Block replacement is convenient, easy to clean and disinfect.



SPECIFICATIONS

Model	BDIB-117
Temperature control range	RT. +5 °C-100 °C
Timer range	1 min-99 h 59 min
Temp. control accuracy	±0.5 °C
Heating time (25 °C to 100 °C)	≤15 min (@20 °C~30 °C)
Certificate	CE
Dimension (LxWxH) mm	225x160x197 mm
Weight	3 kg
Alt Name	Dry Bath Incubator

ACCESSORIES FOR PURCHASE

No	Name	Description
1	Block	96x0.2 mL tube
2	Block	29x0.5 mL+ 32x1.5 mL tube
3	Block	61x0.5 mL tube
4	Block	41x1.5 mL tube
5	Block	41x2.0 mL tube
6	Block	61xΦ6 mm tube
7	Block	28xΦ12 mm tube
8	Block	28xΦ13 mm tube
9	Block	21xΦ15 mm tube
10	Block	21xΦ16 mm tube
11	Block	17xΦ18 mm tube
12	Block	14xΦ20 mm tube
13	Block	9xΦ25 mm tube
14	Block	4xΦ32 mm tube



1



2



3



4



5



6



7



8



9



10



11



12



13



14



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com