



## MINI DRY BATH INCUBATOR BDMI-106

# MINI DRY BATH INCUBATOR BDMI-106

## MINI DRY BATH

A digital Roller Mixer provides a gentle but MINI DRY BATH is a microprocessor -controlled product. It is widely used in preservation and reactions of samples, DNA amplification and initial denaturation of electrophoresis and so forth. It fits almost anywhere and can even be used in cars ,boats or wherever a 24volt power source is available.



Portable;Optional 24V adapter  
LCD display ,it s easy to setup and use  
Compact, fits in the palm of your hand  
Room temperature to 37C in two minutes  
Clear cover ensures temperature uniformity  
Different kinds of aluminum blocks are available

## SPECIFICATIONS

Model	BDMI-106
Temp. Range	+5°C - 100°C
Temp. Accuracy	≤ ±0.3°C
Temp. Uniformity	≤ ±0.3°C
Dimensions (LxWxH, mm)	145 x 113 x 110
Weight	0.5 kg
Electrical	AC 100-240V, 50-60 Hz, 40 W
MOQ-5PCS	
Alt Name	Mini Dry Bath

## ACCESSORIES FOR PURCHASE

No	Name	Specification
1	microtube block	2.0 ml x 15 (MOQ-5PCS)
2	microtube block	1.5 ml x 15(MOQ-5PCS)
3	microtube block	0.5 ml x 24 (MOQ-5PCS)
4	microtube block	0.2 ml x 32 (MOQ-5PCS)



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)