



EXPLOSIVE DRUG DETECTOR BDRD-201

EXPLOSIVE DRUG DETECTOR BDRD-201

EXPLOSIVES AND DRUG DETECTOR

Trace Drug Detector has ion mobility spectrometry (IMS) technology which helps to quickly and accurately detect and determine the types of trace explosives and drugs.



- 2.8 inch TFT touch screen in color display.
- Fast identify speed and accurate result.
- Fashionable design and light weight.
- Top advanced PIMS technology.
- Battery backup time: more than 6 hours.
- High sensitivity: can detect at least 100 Nanogram powder, including fireworks and civil homemade explosives.
- High speed identification.
- Advanced migration tubes ensure high resolution.
- Automatically cleaning systems can prevent pollution.
- Renewable gas purification system.
- Data open to users can add new explosives.

SPECIFICATIONS

Model	BDRD-201
Explosives identify	A variety of military, civilian and homemade explosives: Black Powder, AN, TNT, DNT, Tetra, PETN, Gun Power, NG, RDX, fireworks, C4 etc.
Narcotic detection (option)	Drug, cocaine, opiates (Heroin and Morphine), cannabis (Marijuana and Hashish), and amphetamine-type stimulants (amphetamine, ecstasy & Methamphetamine)
Sampling	Dipstick tracking trace particles
Databases	Open to user, can add new type explosive
Alarm	Audio and visual alert
Sensitivity	Mix: Nanogram explosives, even pictogram sulfur
Sensitivity limit	100 Nanogram TNT
Analysis time	2-10 seconds
Warm-up time	Within 15 minutes
False alarm rate	≤ 1%
Built-in clearing system	Clear the system within 10 seconds
Display	2.8 inch TFT color touch screen
Rechargeable battery standby time	2-4 hours
Operating temperature	-20°C ~ 55°C
Dimensions (LxWxH)	450x160x200 mm
Weight	3.5 kg
Power	AC187-240V 50/60Hz
Alt Name	Explosives And Drug Detector

FEATURES

- 2.8 inch TFT touch screen in color display.
- Fast identify speed and accurate result.
- Fashionable design and light weight.
- Top advanced PIMS technology.
- Battery backup time: more than 6 hours.
- High sensitivity: can detect at least 100 Nanogram powder, including fireworks and civil homemade explosives.
- High speed identification.
- Advanced migration tubes ensure high resolution.

Automatically cleaning systems can prevent pollution.
Renewable gas purification system.
Data open to users can add new explosives.
Dipstick can be reused.
Unlimited storage, USB port to output the data.

APPLICATIONS

Widely used in subway, railway, public security, border police, government Units, postal logistics, other key security zones, the human body clothing, luggage, drugs left on the surface of the goods, explosives, other prohibited items, suspicious objects can be quickly detected.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com