



## BENCHTOP AUTOCLAVE BEM-20202

# BENCHTOP AUTOCLAVE BEM-20202

## TABLE TOP AUTOCLAVE CLASS N SERIES

Class N autoclave discharges the cold air in the sterilization chamber by gravity displacement, saturates the internal steam, and reaches the set pressure and temperature. To achieve the purpose of sterilizing the items to be sterilized. Suitable for stomatology, ophthalmology, beauty, operating room and other places. It is suitable for sterilization of glassware, bare metal utensils, and unencapsulated solid instruments.



High quality stainless steel sterilization chamber and sealing cover.  
Microcomputer control for automatic circulation, with 121°C and 134°C freely settable.  
LED display with touch buttons.  
Equipped with standard test interface.  
European Class N standards.  
Automatic protection functions: over-temperature protection, over-pressure protection, low water level protection to prevent dry burning.  
Buzzer reminder at the end of sterilization and automatic shutdown.  
Automatically discharge cold air. Automatically discharge steam at the end of sterilization.  
Built-in steam generating device to quickly provide saturated steam.

## SPECIFICATIONS

|                           |                                             |
|---------------------------|---------------------------------------------|
| Model                     | BEM-20202                                   |
| Capacity                  | 23L                                         |
| Sterilization Temperature | 121/134°C                                   |
| Working Pressure          | 0.22 MPa                                    |
| Features                  |                                             |
| Door                      | Manual                                      |
| Display                   | LED                                         |
| Max Designed Pressure     | 0.25MPa                                     |
| Max Design Temp.          | 139°C                                       |
| Temp. Precision           | 0.1°C                                       |
| Chamber Material          | SUS304                                      |
| Chamber Size              | φ250*452mm                                  |
| Altitude                  | ≤2500m                                      |
| Rack                      | 1 rack with 2 trays                         |
| Standard Accessories      | Protective gloves, draining hose, door seal |
| Consumption               | 2KW                                         |
| Power Supply              |                                             |
| External Size(W*D*H)      | 600*400*350mm                               |
| Net Weight                | 27kg                                        |
| Package Size(W*D*H)       | 640*450*390mm                               |
| Gross Weight              | 29kg                                        |
| Alt Name                  | Table Top Autoclave Class N Series          |

## FEATURES

High quality stainless steel sterilization chamber and sealing cover.  
Microcomputer control for automatic circulation, with 121°C and 134°C freely settable.  
LED display with touch buttons.  
Equipped with standard test interface.  
European Class N standards.  
Automatic protection functions: over-temperature protection, over-pressure protection, low water level protection to prevent dry burning.

Buzzer reminder at the end of sterilization and automatic shutdown.  
Automatically discharge cold air. Automatically discharge steam at the end of sterilization.  
Built-in steam generating device to quickly provide saturated steam.  
Two sterilization procedures at 121°C and 134°C are selectable.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)