



BENCHTOP AUTOCLAVE BEM-21701

BENCHTOP AUTOCLAVE BEM-21701

HOT AIR STERILIZER



The maximum temperature can reach 250 degrees Celsius to maintain the disinfection effect.
Precise temperature control design can accurately reach the set disinfecting temperature to ensure that the temperature does not drop.
It has an observable perspective window to observe the situation of the disinfected items at any time.
Timing system, hands free, more convenient to use.

SPECIFICATIONS

Model	BEM-21701
Capacity	23L
Sterilization Temperature	5~40°C
Working Pressure	
Features	
Circulation Mode	Forced convection
Temp. Control Way	PT100
Temp. Control Range	RT+10~250°C
Temp. Fluctuation	±3°C
Temp. Uniformity	±5°C
Ambient Temp.	5~40°C
Chamber Material	201 stainless steel
Chamber Size(W*D*H)mm	256*263*314
Altitude	
Shelf(mm)	240*240mm, 2pcs, can be optional
Standard Accessories	Shelf*2
Optional Accessories	Shelf
Consumption	600W
Power Supply	AC220V, 50/60Hz(Standard); AC110V, 50/60Hz(Optional, external transformer)
External Size(W*D*H)	583*393*486mm
Package Size(W*D*H)	675*522*641mm
Net Weight	29.5kg
Gross Weight	40kg
Alt Name	Hot Air Sterilizer

APPLICATIONS

It is used to sterilize items not affected by high temperature, such as glassware, porcelain, gelatin sponge, liquid paraffin, vaseline, ointment, various powders and metal surgical instruments.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com