



GLASS BEAD STERILIZER BEM-21801

GLASS BEAD STERILIZER BEM-21801

GLASS BEAD STERILIZER



Digital control, up to 300°C.

Sterilizes forceps, scissors, needles&other small surgical instruments(generally within 10~15s).

Safe to use: no gases, fumes, open flames or harsh liquids.

Stainless steel inner container, built-in highly precise temperature controller and over-temperature protection.

Includes one bag of glass beads.

SPECIFICATIONS

Model	BEM-21801
Glass Beads Capacity	150g
Temp. Range	100~300°C
Temp. Accuracy	±5°C
Temp. Precision	1°C
Timing Range	1min~99h59min/∞
Heating Time	≤25min(RT. to 300°C)
Material of Lid	Stainless steel
Chamber Size	φ40*80mm
Standard Accessories	Glass beads
Max Consumption	120W
Power Supply	AC110/220V±10%, 50/60Hz
External Size(W*D*H)	176*135*189.5mm
Net Weight	2.5kg
Package Size(W*D*H)	330*300*300mm
Gross Weight	3.8kg
Alt Name	Glass Bead Sterilizer

APPLICATIONS

Glass Bead Sterilizer is ideal for the sterilization of surgical instruments(forceps, scissors, needles, etc.).When the chamber is filled with the included glass beads, high temperature(up to 300°C) can be used to eliminate the presence of bacteria, spores and other microorganisms. Without using potentially dangerous gases, flames or chemicals,it ensures that the sterilized instruments will be free of pathogens and microbial contaminants.The small footprint saves bench space, the sterilizer is safe to be used even in laminar flow hoods.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com