



HYDROGEN PEROXIDE STERILIZER BEM-22801

HYDROGEN PEROXIDE STERILIZER BEM-22801

HYDROGEN PEROXIDE STERILIZER



Efficient broad spectrum sterilization, can kill white staphylococcus, natural bacteria in the air and so on.
Micron atomized particles form Brownian motion, uniform diffusion.
Delay start + wireless remote control one-button to start, safe and convenient.

SPECIFICATIONS

Model	BEM-22801
Display	4.3-inch touch screen
Concentration Range	7.5% (W/W) H ₂ O ₂ disinfectant, 0-188.4 mg/m ³
Atomized Particle Size	1 µm on average
Volume of Spray	15-18 ml/min
Mist Discharge Mode	One-way discharge mist
Working Time	60 min
Stock Solution Volume	1 L
Natural Bacteria Extinction Rate	>90%, wooden countertop extinction log >1
Sterilization Rate of White Staphylococcus Albus	>99.99%
Consumption	1300 W
Power Supply	AC 220 V, 50/60 Hz (Standard); AC 110 V, 50/60 Hz (Optional)
External Size (WxDxH)	452 x 256 x 346 mm
Net Weight	10.5 kg
Package Size (WxDxH)	600 x 326 x 386 mm
Gross Weight	13 kg
Alt Name	Hydrogen Peroxide Sterilizer

APPLICATIONS

It is suitable for biosafety laboratories, infection wards, schools, hospitals, food enterprises, scientific research institutes and other places that require regular disinfection of air and hard surfaces.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com