



HYDROGEN PEROXIDE STERILIZER BEM-22901

HYDROGEN PEROXIDE STERILIZER BEM-22901

H₂O₂ NEBULIZER STERILIZER



It can control the spraying of the disinfectant in the form of nano-micron droplets, reduce the condensed liquid, and ensure the penetration and disinfection effect on complicated and difficult-to-reach places.

Full-color touch screen, easy and clear operation, high user experience comfort. It can achieve a killing effect of 6 logs of *Bacillus subtilis* black mutated spores. It is less corrosive, non-toxic, and has a wide range of applications, which can replace traditional space disinfection methods.

Flexible to use and easy to operate, small in size and easy to move, suitable for various units that have requirements for microbial limits. It has an automatic calculation function. Just input the size of the room to be disinfected, and the required amount of disinfectant can be automatically calculated.

SPECIFICATIONS

Model	BEM-22901
Liquid Capacity	4L
Suitable Space	100 m ³
Consumption	1100 W
Power Supply	AC220V, 50/60Hz (Standard); AC110V, 50/60Hz (Optional) (External transformer)
External Size (WxDxH)	530 x 530 x 950 mm
Net Weight	48 kg
Package Size (WxDxH)	620 x 620 x 960 mm
Gross Weight	55 kg
Alt Name	H ₂ O ₂ Nebulizer Sterilizer

FEATURES

It can control the spraying of the disinfectant in the form of nano-micron droplets, reduce the condensed liquid, and ensure the penetration and disinfection effect on complicated and difficult-to-reach places.

Full-color touch screen, easy and clear operation, high user experience comfort.

It can achieve a killing effect of 6 logs of *Bacillus subtilis* black mutated spores.

It is less corrosive, non-toxic, and has a wide range of applications, which can replace traditional space disinfection methods.

Flexible to use and easy to operate, small in size and easy to move, suitable for various units that have requirements for microbial limits. It has an automatic calculation function. Just input the size of the room to be disinfected, and the required amount of disinfectant can be automatically calculated.

With the function of timing appointments, this product can make an appointment in advance. After the timing appointment is made, the equipment can automatically start disinfection.

APPLICATIONS

It is suitable for biosafety laboratories, infection wards, schools, hospitals, food enterprises, scientific research institutes and other places that require regular disinfection of air and hard surfaces.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com