



VERTICAL LAMINAR AIRFLOW BEM1BE4 (BLVR-204)

VERTICAL LAMINAR AIRFLOW BEM1BE4

VERTICAL LAMINAR FLOW CABINET



SPECIFICATIONS

Model	BEM1BE4
Old Model	BLVR-204
External Size (W*D*H)	1460*620*1850mm
Internal Size W*D*H)	1335*530*650mm
Work Surface Height	750mm
Display	LCD Display
Airflow Velocity	Average of 0.3~0.45m/s, air speed adjustable
Material	Main Body: Cold-rolled steel with anti-bacteria power coatingWork Table: 304 stainless steel
Pre-Filter	Polyester fiber, washable
HEPA Filter	99.995% efficiency at 0.3μm
Noise	≤65dB
Front Window	Manual,5mm toughened glass, anti-UV
Max Opening	400mm
Illuminating Lamp	LED lamp16W*1
UV Lamp	40W*1Emission of 253.7 nanometers
Consumption	250W
Standard Accessory	LED lampUV lamp*2Base stand
Optional Accessory	Electric height adjustable base stand
Caster	Universal wheel with leveling feet
Power Supply	AC220V±10%, 50/60Hz; 110V±10%, 60Hz
Gross Weight	224kg
Package Size (WD'H)	1635*920*1355mm
Alt Name	Vertical Laminar Flow Cabinet



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com