



## PASS BOX BEM1CL3 (BAPB-103)

## PASS BOX BEM1CL3

PASS BOX



## SPECIFICATIONS

Model	BEM1CL3
Old Model	BAPB-103
Type	Air Shower Pass Box
External Size(W*D*H)mm	860*770*1200
Internal Size(W*D*H)mm	700*700*700
Material	304Stainless steel
Interlock	Electronical
UV Lamp	15W*1
HEPA Filter	Efficiency≥99.995% at 0.3μm
Clean Air Speed	0.36~0.54m/s.
Power Supply	220V, 50Hz; 110V, 50/60Hz
Gross Weight (Kg)	120
Packing Size(W*D*H) mm	980*900*1450
Alt Name	Pass Box



## APPLICATIONS

Pass Box is used to transfer goods between rooms of different clean levels of clean room. It can reduce times of opening door, prevent the air convection between different rooms and minimize the pollution extent.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)