



PASS BOX BEM1CL4 (BSPB-101)

PASS BOX BEM1CL4

PASS BOX



SPECIFICATIONS

Model	BEM1CL4
Old Model	BSPB-101
Type	Pass Box
External Size(W*D*H)mm	660*570*580
Internal Size(W*D*H)mm	500*500*500
Material	304Stainless steel
Interlock	Electronical
UV Lamp	8W*1
HEPA Filter	/
Clean Air Speed	/
Power Supply	220V, 50Hz; 110V, 50/60Hz
Gross Weight (Kg)	55
Packing Size(W*D*H) mm	780*700*770
Alt Name	Pass Box



APPLICATIONS

Pass Box is used to transfer goods between rooms of different clean levels of clean room. It can reduce times of opening door, prevent the air convection between different rooms and minimize the pollution extent.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com