



## DOUBLE ADJUSTABLE FREQUENCIES ULTRASONIC CLEANER BEM1CR4

# DOUBLE ADJUSTABLE FREQUENCIES ULTRASONIC CLEANER

## BEM1CR4

### DOUBLE ADJUSTABLE FREQUENCIES ULTRASONIC CLEANER

The Double Adjustable Frequencies Type device is a versatile equipment commonly used in applications requiring precise frequency modulation for various industrial, medical, and research purposes.



1. Over-voltage and over-current protection.
2. Standard: 1pc of basket, 1pc of lid.
3. Low-level water protection.
4. Fully microprocessor controlled and programmable.
5. Digital display, easy setting for time and temperature.
6. Time and temperature real-time display.

## SPECIFICATIONS

|                                       |                                                  |
|---------------------------------------|--------------------------------------------------|
| Model                                 | BEM1CR4                                          |
| Capacity (L)                          | 27                                               |
| Adjustable Ultrasonic Frequency (KHz) | 45/80                                            |
| Time Range / Heating Temp             | 1~480min / 0~80°C                                |
| Ultrasonic Power (W)                  | 700                                              |
| Heating Power (W)                     | 800                                              |
| Tank Size (L*W*H mm)                  | 500*300*180                                      |
| Packing Size (L*W*H mm)               | 1100*690*800                                     |
| Gross Weight (kg)                     | 104                                              |
| Standard Power Supply                 | AC220V±10%, 50/60Hz                              |
| Alt Name                              | Double Adjustable Frequencies Ultrasonic Cleaner |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)