



CLASS I BIOSAFETY CABINET BEM1V2 (BCBS-102)

CLASS I BIOSAFETY CABINET BEM1V2

CLASS I BIOLOGICAL SAFETY CABINET



1. UV lamp for sterilization.
2. Soft touch control panel, LCD display.
3. Energy-saving, high efficiency, low noise.
4. With memory function in case of power-failure.
5. One piece 304 stainless steel work table, easy for cleaning.
6. Stable air flow system; provide firm protection to people and environment.

SPECIFICATIONS

Model	BEM1V2
Old Model	BCBS-102
External Size(W*D*H)	700*550*900mm
Internal Size(W*D*H)	680*540*500mm
HEPA Filter	99.995% efficiency at 0.3 μ m
Airflow Velocity	0.75~1.0m/s
Noise	\leq 55dB
Illuminating Lamp	8W*1LED lamp
UV Lamp	20W*1Emission of 253.7 nanometers for most efficient decontamination
Power Consumption	160W
Work Surface Height	/
Caster	/
Main Body Construction	PMMA
Work Table	304 stainless steel
Power Supply	AC220V \pm 10%, 50/60Hz; 110V \pm 10%, 60Hz
Standard Accessory	LED lamp, UV lamp*2
Optional Accessory	Electric height adjustable base stand
Gross Weight	68kg
Package Size(W*D*H)	870*720*1170mm
Alt Name	Class I Biological Safety Cabinet



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com