



METALLURGICAL MICROSCOPE BEM2B2 (BMIC-802)

METALLURGICAL MICROSCOPE BEM2B2

METALLURGICAL MICROSCOPE

Metallographic microscope is a high-tech product which combines optical microscope technology, photoelectric conversion technology and computer image processing technology perfectly. It can easily observe metallographic images on a computer, analyze and grade metallographic maps, edit, save and print the pictures.



SPECIFICATIONS

Model	BEM2B2
Old Model	BMIC-802
Viewing Head	Binocular Head Inclined at 45°
Nosepiece	Quadruple Nosepiece
Eyepiece	WF10X
Objective	Achromatic Objective: 10X, 100XLWD Plan Achromatic Objective: 40X
Storage	Double Layer Mechanical Stage 200*180mm, qφ10 well table plate, qφ20 well table plate
Focusing	Coaxial coarse&fine focusing adjustment with rack and pinion mechanism fine focusing scale value 0.002mm
Illumination	Build-in Illumination, halogen lamp 6V/20W
Filter	Blue, Green, Yellow
Power Supply	AC110~220V, 50/60Hz
Package Size(mm)	420*350*520
Gross Weight(kg)	10
Alt Name	Metallurgical Microscope



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com