



BIOLOGICAL MICROSCOPE BEM2BE2 (BMIC-504)

BIOLOGICAL MICROSCOPE BEM2BE2

INVERTED MICROSCOPE



SPECIFICATIONS

| | |
|------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Model | BEM2BE2 |
| Old Model | BMIC-504 |
| Viewing Head | Trinocular viewing head, inclined at 30° |
| Eyepiece | EW10X/22mm |
| Objective | LWD Infinite plan objective 4X,10X,20X, 40X Infinite plan objective 10X,20X, 40X |
| Stage | Plain stage: 160*250mm, glass insert, attachable mechanical stage, X-Y coaxial control, moving range: 120*78mm Auxiliary stage 70*180mm, Terasaki holder, Petri dish holder $\phi 35\text{mm}/\phi 54\text{mm}/\phi 90\text{mm}$, slide glass holder $\phi 54\text{mm}$ |
| Condenser | ELWD Condenser NA 0.3, LWD 72mm (without condenser 150mm) |
| Interpupillary Distance | 50-75mm/48-75mm |
| Focusing | Coaxial coarse& fine adjustment, fine division 0.002mm |
| Phase Contrast Sliding Plate | Center of phase contrast ring plate adjustable 10X, 20X, 40X phase annulus plate (fixed) 10X, 20X, 40X phase annulus plate (adjustable) |
| Center Telescope | Center telescope($\phi 30\text{mm}$) |
| Illumination | LED lamp 3W/Halogen lamp 6V/30W |
| Filter | Blue, green and frosted glass |
| Optional Accessories | Inverted fluorescent attachment, photo attachment, video adapter with C mount |
| Power Supply | AC 110/220V $\pm 10\%$,50/60Hz |
| Package Size | 468*397*665mm |
| Gross Weight | 15kg |
| Alt Name | Inverted Microscope |



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com