



BIOLOGICAL MICROSCOPE BEM2BF1 (BMIC-401)

BIOLOGICAL MICROSCOPE BEM2BF1

FLUORESCENCE BIOLOGICAL MICROSCOPE

The Fluorescence Biological Microscope is an advanced imaging instrument designed for high-precision analysis of biological samples under fluorescence illumination.



SPECIFICATIONS

Model	BEM2BF1
Old Model	BMIC-401
View Head	Binocular head
Fluorescent Straight Illumination	B(Blue), G(Green) exciting light filter system O ordinary light system 100W HBO Super high-pressured spherical mercury lamp
Protection Screen	Screen to resist the UV
Nosepiece	Quadruple nosepiece
Connecting Wire	With 2 connectors
Fluorescent Free Objective	4X,10X, 40X, 100X(Oil)
Immersion Oil	Fluorescent free oil
Power Supply	AC100~240V, 50/60Hz
Package Size	320*270*425mm
Gross Weight	7kg
Alt Name	Fluorescence Biological Microscope



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com