



## STEREO ZOOM MICROSCOPE BEM2BH3 (BMIC-905)

# STEREO ZOOM MICROSCOPE BEM2BH3

## STEREO ZOOM MICROSCOPE

A stereo zoom microscope is a kind of microscope that provides stereo vision and continuous zoom function, mainly used for observing and studying objects. It combines the features of stereo microscopy and zoom microscopy to provide clear stereoscopic images in three-dimensional space, and the magnification can be continuously changed by adjusting the lens assembly.



## SPECIFICATIONS

Model	BEM2BH3
Old Model	BMIC-905
Illumination	Natural light
View Head And Body	Inclined viewing body 360°rotatable
Objective	2X&4X, 360° rotated
Working Distance	100mm
Eyepiece	WF10X/20
Power Supply	AC100~240V, 50/60Hz
Package Size(W*D*H)	400*250*360mm
Gross Weight	8kg
Alt Name	Stereo Zoom Microscope



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)