



## STRONG ACID AND ALKALI STORAGE CABINET BEM2C4 (BCBT-204)

# STRONG ACID AND ALKALI STORAGE CABINET BEM2C4

## STRONG ACID & ALKALI STORAGE CABINET

The Strong Acid & Alkali Storage Cabinet is a specialized safety cabinet designed for the safe storage of corrosive chemicals, including strong acids and alkalis, in laboratory and industrial environments.



1. Padlock design.
2. Adjustable shelves.
3. Anti-corrosive materials for strong acid and alkali.
4. Customizable.

## SPECIFICATIONS

Model	BEM2C4
Old Model	BCBT-204
Capacity	28 gallon(105L)
Material	PP
Adjustable Shelf	2pcs
Door	Double/Manual
External Size(W*D*H)mm	910*600*900
Package Size(W*D*H)mm	995*685*1035
Gross Weight	48kg
Alt Name	Strong Acid & Alkali Storage Cabinet



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)