



STRONG ACID AND ALKALI STORAGE CABINET BEM2C6 (BCBT-206)

STRONG ACID AND ALKALI STORAGE CABINET BEM2C6

STRONG ACID & ALKALI STORAGE CABINET

The Strong Acid & Alkali Storage Cabinet is a specialized safety cabinet designed for the safe storage of corrosive chemicals, including strong acids and alkalis, in laboratory and industrial environments.



1. Padlock design.
2. Adjustable shelves.
3. Anti-corrosive materials for strong acid and alkali.
4. Customizable.

SPECIFICATIONS

Model	BEM2C6
Old Model	BCBT-206
Capacity	45 gallon(170L)
Material	PP
Adjustable Shelf	2pcs
Door	Double/Manual
External Size(W*D*H)mm	1090*460*1650
Package Size(W*D*H)mm	1175*545*1785
Gross Weight	72kg
Alt Name	Strong Acid & Alkali Storage Cabinet



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com