



STRONG ACID AND ALKALI STORAGE CABINET BEM2C7 (BCBT-207)

STRONG ACID AND ALKALI STORAGE CABINET BEM2C7

STRONG ACID & ALKALI STORAGE CABINET

The Strong Acid & Alkali Storage Cabinet is a specialized safety cabinet designed for the safe storage of corrosive chemicals, including strong acids and alkalis, in laboratory and industrial environments.



1. Padlock design.
2. Adjustable shelves.
3. Anti-corrosive materials for strong acid and alkali.
4. Customizable.

SPECIFICATIONS

Model	BEM2C7
Old Model	BCBT-207
Capacity	60 gallon(227L)
Material	PP
Adjustable Shelf	2pcs
Door	Double/Manual
External Size(W*D*H)mm	860*860*1650
Package Size(W*D*H)mm	945*945*1785
Gross Weight	90kg
Alt Name	Strong Acid & Alkali Storage Cabinet



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com