



## PORTABLE REFRIGERATOR BEM2CE3 (BPOR-202)

# PORTABLE REFRIGERATOR BEM2CE3

## PORTABLE REFRIGERATOR

A portable refrigerator is a self-contained storage chest which uses an outside power source to control its internal temperature, also it doesn't rely on ice or freezer packs to keep its contents cool.



- 1.The storage container can keep content cooling up to 24 hours
- 2.Insulation material is PU CFC-free polyurethane
- 3.Made of environmental PP material,food grade, non-toxic
- 4.Storage container is durable for vaccine storage, blood products storage,dry ice storage,frozen specimen storage etc.

## SPECIFICATIONS

Model	BEM2CE3
Old Model	BPOR-202
External Size (WxDxH)	440x300x340 mm
Internal Size (WxDxH)	350x220x250 mm
Capacity	17 L
Temperature Range	2°C-8°C
Storage Time (Ambient 30°C)	24 h
Material (External)	PP
Material (Internal)	PU
Ice Vessels	2 (CMB-1200)
Max Loading	15 kg
Optional Accessory	LED display thermometer
Net Weight	3.2 kg
Gross Weight	6.5 kg
Package Size (WxDxH)	460x330x370 mm
Alt Name	Portable Refrigerator

## APPLICATIONS

Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)