



MICROWAVE DIGESTER BEM2CO1 (BMWD-203)

MICROWAVE DIGESTER BEM2C01

MICROWAVE DIGESTER

The chamber is made of stainless steel, resistant to corrosion, and long-lasting use Anti-explosion and anti-corrosion composite fiber outer vessel for aerospace. Safety bolt with multiple uses rather than a safety membrane and other consumables. Exact tension constrained by piezoelectric gem without cross-tainting issue.



All-steel industrial grade chamber, anti-corrosion and durable use.
Aerospace composite fiber outer vessel , anti-explosion and anti-corrosion.
Multi-functional safety bolt design, instead of safety membrane and other consumables.
Precise pressure controlled by piezoelectric crystal without cross-contamination problem.

SPECIFICATIONS

Model	BEM2C01
Old Model	BMWD-203
Vessel Quantity	6
Vessel Volume	100mL
Vessel material	Aerospace composite fiber for outer Vessel, Modified TFM material for inner Vessel
Temperature Range	0-300 °C
Temperature Control Accuracy	±1 °C
Temperature Sensor	High-precision platinum resistor temperature sensor
Pressure Control range	0-10MPa
Pressure Control Accuracy	±0.01MPa
Pressure Sensor	Piezoelectric crystal pressure sensor
Chamber Exhaust System	High-power anti-corrosion axial fan, exhaust speed: 3.1m ³ /min
Digestion Method	Internal software, up to 20 methods can be stored
Display	LCD display
Maximum Output Power	1000W
Power Supply	AC220V ± 10%, 50Hz
External Size(WxDxH)	450x515x510 mm
Package Size(WxDxH)	Main Body: 710x670x700 mm, Accessory: 490x420x310 mm
Gross Weight	Main Body: 63.5kg, Accessory: 9kg
Alt Name	Microwave Digester

APPLICATIONS

Pharmaceutical Analysis, Agricultural, Environmental Analysis, Biological Sciences



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com