



## PORTABLE REFRIGERATOR BEM2F4 (BPOR-102)

## PORTABLE REFRIGERATOR BEM2F4

### PORTABLE REFRIGERATOR

The Portable Refrigerator is a versatile cooling solution designed for travel, camping, medical, and laboratory use, providing precise and reliable temperature control in a compact, durable design. Ideal for storing vaccines, food, and other temperature-sensitive items, this refrigerator combines portability with energy efficiency to ensure stable cooling in any environment.



## SPECIFICATIONS

|                             |                                                                                                   |
|-----------------------------|---------------------------------------------------------------------------------------------------|
| Model                       | BEM2F4                                                                                            |
| Old Model                   | BPOR-102                                                                                          |
| Capacity                    | 12L                                                                                               |
| Internal&External Materials | PP+PE                                                                                             |
| Insulation Material         | PU                                                                                                |
| Temp. Resolution            | 0.1°C                                                                                             |
| Manufacturing Process       | Blow molding+injection molding, burr-free surface                                                 |
| Cooling Source Material     | Sodium polyacrylate                                                                               |
| Cooling Retention Time      | Maintains an internal temperature of 2~8°C for over 24 hours at an ambient temperature up to 43°C |
| Configuration               | Refrigerator*1, phase-change ice box(400ml)*6, digital thermometer*1, shoulder strap*1, tray*1    |
| Internal Size(W*D*H)        | 325*170*225mm                                                                                     |
| External Size(W*D*H)        | 425*262*305mm                                                                                     |
| Package Size(W*D*H)         | 460*355*360mm                                                                                     |
| Net Weight                  | 3.5kg                                                                                             |
| Gross Weight                | 4kg                                                                                               |
| Alt Name                    | Portable Refrigerator                                                                             |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)