



WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET BEM2J1 (BCBT-101)

WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET BEM2J1

WEAK ACID AND ALKALI CHEMICALS STORAGE CABINET

The Weak Acid and Alkali Chemicals Storage Cabinet is a specialized storage solution designed to safely house chemicals with low corrosivity, such as weak acids and weak alkalis.



SPECIFICATIONS

Model	BEM2J1
Old Model	BCBT-101
Capacity	4 gallon(15L)
Material	Double-layer anti-fire cold-rolled steel,PP tray
Adjustable Shelf	1pc
Door	Single/Manual
External Size(W*D*H)mm	430*430*560
Package Size(W*D*H)mm	505*505*585
Gross Weight	28kg
Alt Name	Weak Acid and Alkali Chemicals Storage Cabinet



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com