



## WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET BEM2J4 (BCBT-104)

# WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET

## BEM2J4

### WEAK ACID AND ALKALI CHEMICALS STORAGE CABINET

The Weak Acid and Alkali Chemicals Storage Cabinet is a specialized storage solution designed to safely house chemicals with low corrosivity, such as weak acids and weak alkalis.



## SPECIFICATIONS

Model	BEM2J4
Old Model	BCBT-104
Capacity	30 gallon(114L)
Material	Double-layer anti-fire cold-rolled steel,PP tray
Adjustable Shelf	1pc
Door	Double/Manual
External Size(W*D*H)mm	1090*460*1120
Package Size(W*D*H)mm	1170*540*1260
Gross Weight	101kg
Alt Name	Weak Acid and Alkali Chemicals Storage Cabinet



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)