



WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET BEM2J5 (BCBT-105)

WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET

BEM2J5

WEAK ACID AND ALKALI CHEMICALS STORAGE CABINET

The Weak Acid and Alkali Chemicals Storage Cabinet is a specialized storage solution designed to safely house chemicals with low corrosivity, such as weak acids and weak alkalis.



SPECIFICATIONS

Model	BEM2J5
Old Model	BCBT-105
Capacity	45 gallon(170L)
Material	Double-layer anti-fire cold-rolled steel,PP tray
Adjustable Shelf	2pcs
Door	Double/Manual
External Size(W*D*H)mm	1090*460*1650
Package Size(W*D*H)mm	1170*540*1790
Gross Weight	138kg
Alt Name	Weak Acid and Alkali Chemicals Storage Cabinet



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com