



## WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET BEM2J7 (BCBT-107)

# WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET

## BEM2J7

### WEAK ACID AND ALKALI CHEMICALS STORAGE CABINET

The Weak Acid and Alkali Chemicals Storage Cabinet is a specialized storage solution designed to safely house chemicals with low corrosivity, such as weak acids and weak alkalis.



## SPECIFICATIONS

Model	BEM2J7
Old Model	BCBT-107
Capacity	90 gallon(340L)
Material	Double-layer anti-fire cold-rolled steel,PP tray
Adjustable Shelf	2pcs
Door	Double/Manual
External Size(W*D*H)mm	1090*860*1650
Package Size(W*D*H)mm	1170*940*1790
Gross Weight	200kg
Alt Name	Weak Acid and Alkali Chemicals Storage Cabinet



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)