




BIOSAFETY TRANSPORT BOX BEM2S3 (BTBT-103)

BIOSAFETY TRANSPORT BOX



Model	BEM2S3
Old Model	BTBT-103
Capacity	12L
Internal & External Materials	PP+PE
Insulation Material	PU
Temp. Resolution	0.1℃
Manufacturing Process	Blow molding + injection molding, burr-free surface
Cooling Source Material	Sodium polyacrylate
Cooling Retention Time	Maintains an internal temperature of 2~8℃ for over 24 hours at an ambient temperature up to 43℃
Sponge Test Tube Holder	Standard: 20 holes, compatible with 10ml test tubes; Optional: 11 holes, compatible with 20ml test tubes
Configuration	Box*1, phase-change ice box(400ml)*6, absorbent cotton*2, shoulder strap*1, tray*1, digital thermometer*1, sponge test tube holder*2, 95KPa biosafety transport tank*2
Internal Size(W*D*H)	325*170*225mm
External Size(W*D*H)	425*262*305mm
Package Size(W*D*H)	460*340*340mm
Net Weight	4kg
Gross Weight	4.8kg
Alt Name	Biosafety Transport Box

It can be used for the transportation of biological products requiring physical insulation, such as un2814, un2900, un3373 biological samples, pathogenic microorganism(virus) species, blood, etc.



Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabsscientific.com | Website: www.biolabsscientific.com