



BIOSAFETY TRANSPORT BOX BEM2S5 (BTBT-105)

BIOSAFETY TRANSPORT BOX BEM2S5

BIOSAFETY TRANSPORT BOX



1. Real time temperature value display
2. Refrigeration and temperature control to ensure the consistency of samples
3. Easy to operate and could be used repeatedly
4. Meet the basic requirements of cold chain transportation and storage of vaccines

SPECIFICATIONS

Model	BEM2S5
Old Model	BTBT-105
Capacity	20L
Internal & External Materials	PP+PE
Insulation Material	PU
Temp. Resolution	0.1°C
Manufacturing Process	Blow molding + injection molding, burr-free surface
Cooling Source Material	Sodium polyacrylate
Cooling Retention Time	Maintains an internal temperature of 2~8°C for over 24 hours at an ambient temperature up to 43°C
Sponge Test Tube Holder	Standard: 20 holes, compatible with 10ml test tubes; Optional: 11 holes, compatible with 20ml test tubes
Configuration	Box*1, phase-change ice box(400ml)*9, shoulder strap*1, absorbent cotton*4, digital thermometer*1, tray*1, divider panel*1, sponge test tube holder*4, 95KPa biosafety transport tank*4
Internal Size(W*D*H)	380*245*230mm
External Size(W*D*H)	450*330*300mm
Package Size(W*D*H)	580*320*380mm
Net Weight	4.8kg
Gross Weight	5.5kg
Alt Name	Biosafety Transport Box

APPLICATIONS

It can be used for the transportation of biological products requiring physical insulation, such as un2814, un2900, un3373 biological samples, pathogenic microorganism(virus) species, blood, etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com