



## BIOLOGICAL MICROSCOPE BEM2Z4 (BMIC-704)

# BIOLOGICAL MICROSCOPE BEM2Z4

## ECONOMICAL BIOLOGICAL MICROSCOPE

The Economical Biological Microscope is an affordable and reliable solution for educational, laboratory, and clinical applications. Designed to meet the needs of students, researchers, and professionals, this microscope provides clear and detailed imaging of biological specimens at an exceptional value.



## SPECIFICATIONS

Model	BEM2Z4
Old Model	BMIC-704
Viewing Head	Dual viewing head, one tube 30° inclined and another tube vertical, 360° rotatable
Achromatic Objectives	4X, 10X, 40X(Spring), 100X(Spring, Oil)
Eyepiece	Wide field eyepiece: WF10X(WF16X optional)
Stage	Double layer mechanical stage size 115*125mm
Focusing	Coaxial coarse and fine adjustment, focusing range 30mm, focusing interval 0.002mm
Condenser	Abbe NA=1.25 with Iris diaphragm&filter
Illumination	LED(1W) lamp 220V or 110V(LED 3W lamp optional)
Power Supply	AC110/220V±10%, 50/60Hz
Package Size	450*360*240mm
Gross Weight	5kg
Alt Name	Economical Biological Microscope



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)