



## MAGNETIC STIRRER BER-4202

# MAGNETIC STIRRER BER-4202

## MINI MAGNETIC STIRRER



1. Shell is all plastic material (ABS), lightweight and durable.
2. Mixing smoothly and fast, stirring speed stepless adjustable

## SPECIFICATIONS

Model	BER-4202
Mode	Without light
Function	
Max. Volume (ml)	1000
Speed (r.p.m)	100 ~ 500 or 100 ~ 2000
Working Time	Continuous
Minimum Resolution (R.P.M.)	100
Structure	
Exterior	Plastic
Driving Mode	Motor
Motor Power	0.01kW
Controller	
Continuous Operation	Yes
Setting Mode	Scale
Speed Display Mode	Scale
Timer	No
Operation Function	No
Speed & Temp. Sensor	No
Specification	
Exterior Size (W*L*H)(mm)	φ 138*60
Packing Size (W*L*H)(mm)	210*170*88
Plate Size (mm)	φ 128
Power Supply (50/60Hz) Current Rating	100V~240V
NW/GW (kg)	0.6/0.7
LED Light	Without LED light
Alt Name	Mini Magnetic Stirrer





**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)