




MULTISTATION STIRRER BER-4301

MULTISTATION STIRRER BER-4301

FOUR ROWS MAGNETIC STIRRER



1. 4 unit stainless hotplate can work together or separately.
2. Stirring speed is stepless adjusted.
3. Magnetic mixing equipment.

SPECIFICATIONS

Model	BER-4301
Function	
Max. Volume (ml)	1000ml*4
Speed (r.p.m)	100 ~ 2000
Working Time	Continuous
Minimum Resolution	100
Structure	
Exterior	Cold plate spraying, stainless steel plate
Driving Mode	Motor
Motor Power	0.01kW
Controller	
Continuous Operation	Yes
Setting Mode	Scale
Speed Display Mode	Scale
Timer	No
Operation Function	No
Specification	
Exterior Size (W*L*H)(mm)	300*270*80
Packing Size (W*L*H)(mm)	340*310*110
Plate Size (mm)	300*250
Plate Load	5KG
Power supply (50/60Hz) Current Rating	110 or 220V
NW/GW (kg)	3/3.6
Alt Name	Four Rows Magnetic stirrer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com