



## HEATING MANTLE BER-4505

# HEATING MANTLE BER-4505

## MAGNETIC STIRRING HEATING MANTLES



1. Heating and stirring can be used simultaneously, stepless speed adjustable.
2. Up to 450° T general heating equipment.

## SPECIFICATIONS

Model	BER-4505
Function	
Capacity	1000ml
Heating Mode	Surface heat-conduction
Surface Max.Temp.(°C)	450°C
Speed Range	0 - 1600rpm
Structure	
Heating Element	Nichrome wire
Stir Bar	PTFE
Insulation Layer	High quality aluminum silicate cotton
Shell	Cold roller steel plate, surface spraying
Motor Power	0.01kW
Heating Power (kW)	0.35
Controller	
Working Mode	Continuous
Specification	
Exterior Size (mm)	Φ280x220
Packing Size (mm)	290x290x300
Voltage(V)	100-120V or 200-240V
N.W. (Kg)	5.5
G.W. (Kg)	6
Alt Name	Magnetic Stirring Heating Mantles



## MORE INFO

We can produce different size heating mantle as customer' s requirement. For example: 4000ml,6000ml,50 liters, 100 liters and etc.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)