



## HOTPLATE BER-5003

# HOTPLATE BER-5003


NEW HOT PLATE, GRAPHITE HOT PLATE



1. hot plate or Graphite hot plate, uniform surface temperature.
2. up fast and uniform.
3. Digital hot plate: digital display indicating temperature of hot plate clearly.

## SPECIFICATIONS

|                            |                                      |
|----------------------------|--------------------------------------|
| Model                      | BER-5003                             |
| Function                   |                                      |
| Max Temp. (°C)             | 400°C                                |
| Structure                  |                                      |
| Outer Shell                | Iron plate welding, surface spraying |
| Hot plate                  | Graphite Plate                       |
| Heating Power (W)          | 1800                                 |
| Setting Method             | Digital                              |
| Timer                      | No / Digital                         |
| Temperature Sensor         | K model                              |
| Temp. Accuracy             | ±1%                                  |
| Plate Temp. Uniformity     | ±10%                                 |
| Specification              |                                      |
| Exterior Size (W*L*H) (mm) | 340*370*220                          |
| Packing Size (W*L*H) (mm)  | 550*410*270                          |
| Plate Size (mm)            | 300*300                              |
| Plate Load                 | 20KG                                 |
| Power Supply (50/60Hz)     | 220V                                 |
| NW/GW (kg)                 | 12/13.5                              |
| Alt Name                   | Graphite Hot Plate                   |





**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)