



## HOTPLATE BER-5006

# HOTPLATE BER-5006

NEW HOT PLATE, GRAPHITE HOT PLATE



1. hot plate or Graphite hot plate, uniform surface temperature.
2. up fast and uniform.
3. Digital hot plate: digital display indicating temperature of hot plate clearly.

## SPECIFICATIONS

Model	BER-5006
Function	
Max Temp. (°C)	400°C
Structure	
Outer Shell	Iron plate welding, surface spraying
Hot plate	Graphite Plate
Heating Power (W)	3000
Setting Method	Digital
Timer	No / Digital
Temperature Sensor	K model
Temp. Accuracy	±1%
Plate Temp. Uniformity	±10%
Specification	
Exterior Size (W*L*H) (mm)	420*500*220
Packing Size (W*L*H) (mm)	600*560*270
Plate Size (mm)	450*350
Plate Load	30KG
Power Supply (50/60Hz)	220V
NW/GW (kg)	19/20
Alt Name	Graphite Hot Plate





**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)