



HOTPLATE BER-5008

HOTPLATE BER-5008


NEW HOT PLATE, GRAPHITE HOT PLATE



1. hot plate or Graphite hot plate, uniform surface temperature.
2. up fast and uniform.
3. Digital hot plate: digital display indicating temperature of hot plate clearly.

SPECIFICATIONS

Model	BER-5008
Function	
Max Temp. (°C)	350°C
Structure	
Outer Shell	Iron plate welding, surface spraying
Hot plate	Aluminum tensile
Heating Power (W)	3500
Setting Method	Digital
Timer	No / Digital
Temperature Sensor	K model
Temp. Accuracy	±1%
Plate Temp. Uniformity	±10%
Specification	
Exterior Size (W*L*H) (mm)	470*650*220
Packing Size (W*L*H) (mm)	660*620*250
Plate Size (mm)	600*400
Plate Load	40KG
Power Supply (50/60Hz)	220V
NW/GW (kg)	15/17
Plug	without plug.
Alt Name	Hot Plate





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com