



## WATER BATH BER-5701

# WATER BATH BER-5701

## BIG SIZE WATER BATHS



1. The enclosure of the product is formed and machined by using high-quality steel plate. Static electric spraying process is adopted on the surface, which is sturdy and durable. The inner container is finished by the stainless steel stretching.
2. The liner and upper cover are made of high-quality stainless steel plate, featuring strong corrosion resistance.
3. -shaped heating pipe is adopted for direct heating in water. The temperature rise is quick and the thermal loss is small.
4. Single-row digital display or intelligent temperature controller boasts simple operation and favorable application effect.

## SPECIFICATIONS

Model	BER-5701
Supply voltage	220V/60HZ
Power(W)	2.4kw
Temp-motion	+0.5°C
Temp-range	RT+5 ~ 100°C
Senility of Temp Control	≤± 1
Display error	≤± 2.5
Chamber Size(mm)	620*500*200mm
Products Size(mm)	681*560*301mm
Packing Size(mm)	595*388*445mm
N.W.(Kg)	15KG
G.W.(Kg)	20KG
Alt Name	Big Size Water Bath



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)