



## ULTRASONIC CLEANER BER-5901

# ULTRASONIC CLEANER BER-5901

ULTRASONIC CLEANER DIGITAL SERIES, WITH TIMER



1. just tap water, or industrial alcohol and solvent cleaner for more higher cleaning requirement.
2. timer, clear digital LCD display.
3. Power supply can be either AC 100 120V or AC 200 240V.
4. control chip microcontroller. flexible circuit boards control, more secure & stable.
5. Stainless steel made, resistance to wear.
6. Warranty time: 12 months.

## SPECIFICATIONS

Model	BER-5901
Frequency	42KHz
Transducer (PCS)	1
Capacity	0.8
Ultrasonic Power(W)	50
Timer	0-30 min
Tank Size (L*W*H MM)	150*85*65
Unit Size (L*W*H MM)	180*110*150
G.W/N.W	1.3/1
Alt Name	Ultrasonic Cleaner



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)