



PORTABLE STERILIZER BER-6101

PORTABLE STERILIZER BER-6101

PORTABLE PRESSURE STEAM STERILIZER



1. Fully stainless steel structure.
2. Double scale indication pressure gauge.
3. Electric heated.
4. Easy to operate, safe and reliable.

SPECIFICATIONS

Model	BER-6101
Capacity	18L ϕ 280*260mm
Working Pressure	0.14-0.16Mpa
Working Temperature	126°C
Max.Safety Pressure	0.165Mpa
Heat Average	$\leq \pm 0.1^\circ\text{C}$
Power	2KW/AC200V.50Hz
Packing Size	410*410*520mm
G.W/N.W	18kg/16kg
Alt Name	Portable Pressure Steam Sterilizer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com