



PORTABLE STERILIZER BER-6102

PORTABLE STERILIZER BER-6102

PORTABLE PRESSURE STEAM STERILIZER



1. Fully stainless steel structure.
2. Double scale indication pressure gauge.
3. Electric heated.
4. Easy to operate, safe and reliable.

SPECIFICATIONS

Model	BER-6102
Capacity	24LΦ280*390mm
Working Pressure	0.14-0.16Mpa
Working Temperature	126°C
Max.Safety Pressure	0.165Mpa
Heat Average	± 0.1°C
Power	2KW/AC200V.50Hz
Packing Size	410*410*630mm
G.W/N.W	19kg/17kg
Alt Name	Portable Pressure Steam Sterilizer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com