



## PORTABLE STERILIZER BER-6302

# PORTABLE STERILIZER BER-6302

PORTABLE PRESSURE STEAM STERILIZER



1. Fully stainless steel structure.
2. Hand wheel type of quick-open door structure.
3. Indicator light indicates working state.
4. Double scale indication pressure gauge.
5. Time scope: 0~60 min.
6. Automatically shut off with beep reminding after sterilization.
7. heated.
8. Easy to operate, safe and reliable.

## SPECIFICATIONS

Model	BER-6302
Capacity	24L
Working Pressure	0.14-0.16Mpa
Working Temperature	126°C
Temp. range	126°C
Max.Safety Pressure	0.165Mpa
Timer range	0-60min
Heat average	$\leq \pm 1^{\circ}\text{C}$
Power	AC220V.50Hz/2Kw
Packing Size(mm)	410*410*550
G.W/N.W	17/15KG
Alt Name	Portable Pressure Steam Sterilizer



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)