



## ACCESSORIES BER-6601

# ACCESSORIES BER-6601

BASIC POWER SUPPLY



1. Molding shell, touch keys, microprocessor intelligent control.
2. With automatic memory function.
3. With standard, timing operation function.
4. Automatic detection of no-load, over-load.

## SPECIFICATIONS

Model	BER-6601
Type of Output	Constant-Voltage, Current
Output range	5~300V, 1~300mA, Maximum 90W
Increment	1V, 1mA
Timer Range	1min~9hr, 59min
Display	LED
Output jacks	Two sets of output jacks
Alt Name	Basic Power Supply



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)