



ACCESSORIES BER-6701

ACCESSORIES BER-6701

BASIC POWER SUPPLY



1. Molding shell,touch keys,microprocessor intelligent control.
2. With automatic memory function.
3. With standard,timing operation function.
4. Automatic detection of no-load,over-load.

SPECIFICATIONS

Model	BER-6701
Type of Output	Constant-Voltage, Current
Output range	5~300V,1~300mA,Maximum 90W
Increment	1V,1mA
Timer Range	1min~9hr,59min
Display	LED
Output jacks	Two sets of output jacks
Alt Name	Basic Power Supply



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com