



ACCESSORIES BER-6801

ACCESSORIES BER-6801

UNIVERSAL POWER SUPPLY



1. Molding shell, touch keys, microprocessor intelligent control.
2. With automatic memory function.
3. With standard, timing operation function.
4. Automatic detection of no-load, over-load.

SPECIFICATIONS

Model	BER-6801
Type of Output	Constant-Voltage, Current or Power
Output range	5~300V, 1~400mA 1~120W
Increment	1V, 1mA
Timer Range	1min~99hr, 59min
Display	LCD display with backlight
Output jacks	Four sets of output jacks
Alt Name	Universal Power Supply



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com