



ACCESSORIES BER-6901

ACCESSORIES BER-6901

UNIVERSAL POWER SUPPLY



1. Molding shell, touch keys, dual core microprocessor intelligent control
2. can store 10 electrophoresis methods; With automatic memory function.
3. With constant voltage, constant current, constant power, misoperation, fault intelligent prompt functions.
4. Automatic detection of no-load, over-load, short circuit, rapid resistance change, ground leak and system overheating.

SPECIFICATIONS

Model	BER-6901
Type of Output	Constant-Voltage, Current or Power
Output range	10~600V,1~500mA1~300W
Increment	1V,1mA
Timer Range	1min~99hr,59min
Display	LCD display with backlight
Output jacks	Four sets of output jacks
Alt Name	Universal Power Supply



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com