



HORIZONTAL GEL ELECTROPHORESIS BER-7001

HORIZONTAL GEL ELECTROPHORESIS BER-7001

HORIZONTAL ELECTROPHORESIS CELL



1. Abundant buffer not only guarantees cooling effect, but also keeps the pH Value stable during the whole experiment process.
2. Unfolded electrodes make repair and replacement easy and quick.
3. With the sample background color design on the tray, convenient to adding the sample.
4. Safety cover open button design, convenient opening the cover.
5. The transparent top lid effectively prevents liquid in the cell from volatilizing and avoids and electricity leakage.
6. Auto-switch-off when the lid is removed.

SPECIFICATIONS

Model	BER-7001
Gel Size (W*L)	48*75mm
Sample volume	3.5,9.1
Buffer Volume	~150ml
Dimension (L*W*H)	250*90*50(mm)
Net weight	0.4kg
Alt Name	Horizontal Electrophoresis Cell



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com