


# HORIZONTAL GEL ELECTROPHORESIS BER-7101

# HORIZONTAL GEL ELECTROPHORESIS BER-7101

## HORIZONTAL ELECTROPHORESIS CELL



1. Abundant buffer not only guarantees cooling effect, but also keeps the pH Value stable during the whole experiment process.
2. Unfolded electrodes make repair and replacement easy and quick.
3. With the sample background color design on the tray, convenient to adding the sample.
4. The transparent top lid effectively prevents liquid in the cell from volatilizing and avoids electricity leakage.
5. Auto-switch-off when the lid is removed.

## SPECIFICATIONS

Model	BER-7101
Gel Size (W*L)	78*100mm
Sample volume	10,15
Buffer Volume	~400ml
Dimension (L*W*H)	365*105*55(mm)
Net weight	0.6kg
Alt Name	Horizontal Electrophoresis Cell



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)