



# LABORATORY LOW SPEED CENTRIFUGE BER-7501

# LABORATORY LOW SPEED CENTRIFUGE BER-7501

DIGITAL TABLETOP LOW SPEED CENTRIFUGE



1. DC brushless motor drive, digital display speed and running time
2. Mechanical door lock
3. Used for fat graft purification, PRP and other cosmetic surgery fields.

## SPECIFICATIONS

Model	BER-7501
Maximum speed	4000r/min
Maximum RCF	2350 x g
Max volume	12 x 15ml
Power supply	AC 220V 50HZ
Timer	1 ~ 99min
Noise	≤55dBA
Exterior Size (mm)	270 x 290 x 230
Packing Size (mm)	350 x 340 x 360
N.W / G.W (KG)	8 / 9
Alt Name	Digital Tabletop Low Speed Centrifuge



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)