



MAGNETIC STIRRER BER1BQ2

MAGNETIC STIRRER BER1BQ2

MINI MAGNETIC STIRRER



1. Shell is all plastic material (ABS), lightweight and durable.
2. Mixing smoothly and fast, stirring speed stepless adjustable
1. Shell is all plastic material (ABS), lightweight and durable.
2. Mixing smoothly and fast, stirring speed stepless adjustable
3. Without LED light

SPECIFICATIONS

Model	BER1BQ2
Mode	Without light
Function	
Max. Volume (ml)	1000
Speed (r.p.m)	100 ~ 500 or 100 ~ 2000
Working Time	Continuous
Minimum Resolution (R.P.M.)	100
Structure	
Exterior	Plastic
Driving Mode	Motor
Motor Power	0.01kW
Controller	
Continuous Operation	Yes
Setting Mode	Scale
Speed Display Mode	Scale
Timer	No
Operation Function	No
Speed & Temp. Sensor	No
Specification	
Exterior Size (W*L*H)(mm)	φ 138*60
Packing Size (W*L*H)(mm)	210*170*88
Plate Size (mm)	φ 128
Power Supply (50/60Hz) Current Rating	100V~240V
NW/GW (kg)	0.6/0.7
Alt Name	Mini Magnetic Stirrer





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com