



MULTISTATION STIRRER BER1BR2

MULTISTATION STIRRER BER1BR2

FOUR ROWS MAGNETIC STIRRER



1. 4 unit stainless hotplate can work together or separately.
2. Stirring speed is stepless adjusted.
3. Magnetic mixing equipment.

SPECIFICATIONS

Model	BER1BR2
Function	
Max. Volume (ml)	3000ml*4
Speed (r.p.m)	100 ~ 2000
Working Time	Continuous
Minimum Resolution	100
Structure	
Exterior	Cold plate spraying, stainless steel plate
Driving Mode	Motor
Motor Power	0.02kW
Controller	
Continuous Operation	Yes
Setting Mode	Scale
Speed Display Mode	Scale
Timer	No
Operation Function	No
Specification	
Exterior Size (W*L*H)(mm)	400*320*120
Packing Size (W*L*H)(mm)	440*360*150
Plate Size (mm)	400*300
Plate Load	20KG
Power supply (50/60Hz) Current Rating	110 or 220V
NW/GW (kg)	6/7
Alt Name	Four Rows Magnetic stirrer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com